

# Guildford FNA Needle

The Guildford FNA Needle is specifically designed to improve the effectiveness of the fine needle aspiration procedure and has been designed in association with the Jarvis Breast Screening Centre in Guildford, UK. The unique needle design enables easy progress through tissue, whilst the long fine point reveals a greater cross-section for increased cell collection. The needle length is made to an accuracy of 0.5mm - a much tighter tolerance than standard needles; this makes it ideal for use with stereotaxic machines.



Code	Product Description	Length	Pack Quantity
8500	Guildford FNA Needle 21g	80 mm	10
8510	Guildford FNA Needle 21g	100 mm	10