




Mediplus Mediloops are specially designed to help identify and retract vessels during surgical procedures. Mediloops are made from Silicone Elastomer and are ideal for temporary ligation, occlusion and identification of vessels and structures.

- Oval shape and elasticity allows gentle control without binding or cutting
- Mediloops remain soft and smooth during long procedures - ensuring quick release without risk of adhesion
- Two dimensions are available - for ease of use for different procedures
- Non-porous, so will not stain or obscure
- Colour coded - for ready identification of arteries, veins, ureters and tendons
- Sterile, non-toxic silicone
- Radiopaque
- Specialised packaging prevents kinking and allows easy removal



7014 - Mediloops
Length 40cm



7018/7028 - Mini Mediloops
Length 40cm

Code	Product Description	Pack Quantity
MEDILOOPS		
7014 Mixed	5 loops - one of each colour per pouch	10 pouches
7014 Yellow	5 loops of yellow per pouch	10 pouches
7014 Green	5 loops of green per pouch	10 pouches
7014 Red	5 loops of red per pouch	10 pouches
7014 White	5 loops of white per pouch	10 pouches
7014 Blue	5 loops of blue per pouch	10 pouches
MINI MEDILOOPS		
7018 Yellow	5 loops of yellow per pouch	10 pouches
7028 Red	2 loops of red per pouch	10 pouches

T7