

Safe-Gel Desiccant Sachets for Solution Control

To protect the healthcare profesional from transmission of diseases such as AIDS and Hepatitis B/C through contact with blood and body fluids.



Prevents Spillage



Prevents Leakage



Prevents Overflow

- **SAFE-GEL** is a unique water-soluble film sachet to gel up fluids (up to 500ml per sachet)
- **SAFE-GEL** is a safe, cost-effective unit dosing system. Sachets can be used before or after procedure, and when used after an operative procedure comparisons can be made of input/output values of fluids in TURPs, TCRES
- **SAFE-GEL** reduces risk of contact with blood or body fluids
- **SAFE-GEL** allows implementation of Universal Precautions procedures
- **SAFE-GEL** can be used in:

Operating Departments - in Suction Receivers

Intensive Care Units - in Suction Receivers and Chest Drain

Wards - in disposable Urinals and Bedpans

Nursing Homes - in Commodes, Vomit Bowls

Anywhere a need exists to control flow or liquids to prevent leakage, spillage or overflow

Code	Product Description	Pack Quantity
7500	Safe-Gel desiccant sachets	200

Note: Safe-Gel should not be used before an operative procedure when suction jars are used in tandem
Sachets will not rupture if solvent or oil-based liquids are present
Sachets will dissolve slowly in the presence of whole blood.
COSHH Data available upon request