

# Cystometry Sets

Code	Product Description	Length	Pack Quantity
<b>Cystometry Sets (Latex Free) - Sterile</b>			
5000	Cystometry Set Comprising 1 each - [Code 5100] 4.5 Fr./Ch. Bladder Pressure Catheter, PVC [Code 5200] 10 Fr./Ch. Bladder Filling Catheter, PVC [Code 5415i] 4.5 Fr./Ch. Abdominal Pressure Catheter, PVC Balloon, with introducer	200cm 40cm 200cm	10
5102	Cystometry Set Comprising 1 each - [Code 5100] 4.5 Fr./Ch. Bladder Pressure Catheter, PVC [Code 5200] 10 Fr./Ch. Bladder Filling Catheter, PVC	200cm 40cm	10
5103	Cystometry Set & Rectal Catheter Set Comprising 1 each - [Code 5100] 4.5 Fr./Ch. Bladder Pressure Catheter, PVC [Code 5200] 10 Fr./Ch. Bladder Filling Catheter, PVC [Code 5415] 4.5 Fr./Ch. Abdominal Pressure Catheter, PVC Balloon	200cm 40cm 200cm	10
5104	Cystometry Set & Rectal Catheter Introducer Set Comprising 1 each - [Code 5716] 6 Fr./Ch. Double Lumen Cystometry Catheter [Code 5415j] 4.5 Fr./Ch. Abdominal Pressure Catheter, PVC Balloon, with introducer	65cm* 200cm	10
5105	Cystometry Set & Rectal Catheter Introducer Set Comprising 1 each - [Code 5702] 8 Fr./Ch. Double Lumen Cystometry Catheter [Code 5415j] 4.5 Fr./Ch. Abdominal Pressure Catheter, PVC Balloon, with introducer	65cm* 200cm	10
5106	Cystometry Set & Rectal Catheter Set Comprising 1 each - [Code 5716] 6 Fr./Ch. Double Lumen Cystometry Catheter [Code 5415] 4.5 Fr./Ch. Abdominal Pressure Catheter, PVC Balloon	65cm* 200cm	10
5108	Cystometry Set & Rectal Catheter Set Comprising 1 each - [Code 5702] 8 Fr./Ch. Double Lumen Cystometry Catheter [Code 5415] 4.5 Fr./Ch. Abdominal Pressure Catheter, PVC Balloon	65cm* 200cm	10

\* NB: Pressure line provides additional 130cm length  
See pages U1 & U2