









Pump Tube Sets

A well designed Pump Tube Set is an important component in the smooth running of urodynamics tests such as filling and voiding cystometry. Mediplus Pump Tube Sets are designed to work with the appropriate peristaltic pump to deliver fluid at a predetermined flow rate. Each set terminates in a rotating male luer connector for secure fitment to urodynamic catheters and offers excellent value for money.

Code	Product Description	Pack Quantity
Pump Tube Set for use with following Urodynamics System (Latex Free)		
5801/1	Lectromed 	25
5806	Laborie 	25
5820	MMS UD2000 	25
5822	MMS Libra, Libraplus & Solar 	25
5830	Lifetech Urojanus 	25
5840	Sedia 	25
5872	Medtronic/Dantec 	25
993	Albyn 	25

Also available:

5001	for Mediwatch/Medtronic/Dantec Urodynamics System, with vented spike	25
9901	for Micromedics urodynamics system	20
5850	Gravity Fill Set, 1.8m, luer lock connector	50
5851	Gravity Fill Set, 1.8m, fir tree connector	50
5180	Damping Tube Set	50
5808	Pump Tube Set with Y-connector for Laborie Urodynamics System	25
D00420	Pump Tube for Ellipse/Helix Urodynamics System	10
D00425	Pump Tube Extension	25
D00501	Tube Set for 3 External Pressure Sensors	10