


SetGuard & Extension Sets

SetGuard is a short connecting tube incorporating a non-return valve, designed to fit between a pump tube set and filling catheter, which ensures the pump tube set remains free from contamination caused by back-flow from the catheter. Supplied sterile for single-use, SetGuard enables pump tube sets to be used more than once *safely and with no risk to patients*.

CystoSetGuard, an innovative product for cystoscopy, is a simple extension set incorporating a one-way valve to prevent backflow into cystoscopy tubing, facilitating safe multiple use of the fluid and administration sets. Designed for ease of use and to make a significant cost-saving on cystoscopy consumables, it also reduces waste and makes environmental savings on urology disposables.

Our extension sets are supplied in a variety of lengths and connectors to enable use in most situations with most appliances.



Code	Product Description	Tubing i.d./o.d.	Length	Pack Quantity
SetGuard & CystoSetGuard (<i>Latex Free</i>)				
5182	SetGuard	5mm / 6.8mm	30cm	25
9806	CystoSetGuard	5mm / 6.8mm	35cm	50
9807	CystoSetGuard - extended	5mm / 6.8mm	57cm	25
Extension Sets (<i>Latex Free</i>)				
9003	PVC Extension Set	1.6mm / 3.17mm	5cm	25
9015	PVC Extension Set	1.6mm / 3.17mm	15cm	50
9025	PVC Extension Set	1.6mm / 3.17mm	25cm	50
9050	PVC Extension Set	1.6mm / 3.17mm	50cm	50
9100	PVC Extension Set	1.6mm / 3.17mm	100cm	50
91250	PVC Extension Set	1.7mm / 3.0mm	125cm	50
91255	PVC Extension Set	1.7mm / 3.0mm	125cm	50
9200	PVC Extension Set	1.6mm / 3.17mm	200cm	50
9102	Drainage Bag Adaptor with funnel	4.8mm / 6.8mm	60cm	10