


Single Lumen Filling & Cystometry

Mediplus bladder pressure catheters are 200cm long, enabling direct connection between the patient's bladder and the pressure transducers, so saving time and expense. Each filling catheter in the range has a diamond-shaped eye to enable easy disengagement of the pressure catheter following insertion.


All catheters are fitted with winged luer connectors for secure and easy attachment to and disconnection from the transducer or filling line.

For greater user convenience we also offer Cystometry Sets comprising bladder pressure and filling catheters or bladder filling, pressure and rectal catheters in one pack.

Bladder Pressure Catheter



Bladder Filling Catheter



Code	Product Description	Length	Pack Quantity
Bladder Pressure Catheters (<i>Latex Free</i>)			
5100N	3.3 Fr./Ch. NYLON	200cm	25
5100	4.5 Fr./Ch. PVC	200cm	25
5110	6 Fr./Ch. PVC	200cm	25
Bladder Filling Catheters (<i>Latex Free</i>)			
5206	6 Fr./Ch. PVC	40cm	10
5208	8 Fr./Ch. PVC	40cm	10
5200	10 Fr./Ch. PVC	40cm	10
5200S	10 Fr./Ch. PVC (Soft)	40cm	10
5212	12 Fr./Ch. PVC	40cm	10