


Suprapubic Cystometry Catheter Set

Utilising the Seldinger catheter introduction system to place the custom designed cystometry catheter offers a number of advantages to both the clinician and the patient:

- Minimum trauma for the patient
- Minimum tissue damage
- Locating the bladder with the 18g hypodermic needle confirms correct location for catheter introduction
- The 70%/30% split lumen configuration ensures a faster filling rate, and better pressure measurement using the smallest charriere size
- Colour-coded pressure lines enable easy identification



Code	Product Description	Length	Pack Quantity
Suprapubic Cystometry Catheter Sets (<i>Latex Free</i>)			
5752	6 Fr./Ch.	16cm	5
5750	8 Fr./Ch.	19cm	5