


Other Products Available


9879 Urology Y-Set with one-way valve

Fluid administration set for delivering fluid into the bladder from two containers. For flexi or ridged cystoscopy.



5603 Flow Control Set

For UDS or GI, which provides accurate control of the amount of fluid flow from 0 to 350 millilitres per minute.




Code	Product Description	Pack Quantity
9879	Urology Y-Set with one-way valve	20
5603	Flow Control Set	20

Other Products Available

5600 Fir Tree Female Luer Connector

The adaptor enables the connection of a foley catheter or a neleton catheter to be converted into a female luer. Useful when giving set or drainage aystem terminates in male luer.



Code	Product Description	Pack Quantity
5600	Fir Tree Female Luer Connector	25
MX491P1	Dual Function Luer Lock Cap for male/female ports, White (Non-PVC) *Minimum order quantity 5 boxes (1 x 5 x 100)	100*