


SiVax Silicone Haematuria Catheters for Prostatectomy

The Mediplus SiVax Haematuria Catheter is made of transparent silicone to enable easy visual inspection, with low profile, uniformly inflating balloon to ensure ease of insertion and removal and bladder comfort while in situ. A radiopaque line allows x-ray confirmation of placement, and a retained Luer cap is fitted (to prevent loss of the cap).

Available in four tip styles, in the most widely used sizes, Si-Vax catheters offer a cost-effective solution to your Haematuria needs following Prostatectomy.



| Code | Product Description | Pack Quantity |
|------|---|---------------|
| 4020 | SiVax Silicone Haematuria Catheter, 3-way, 20Fr. - Straight Round Tip | 10 |
| 4022 | SiVax Silicone Haematuria Catheter, 3-way, 22Fr. - Straight Round Tip | 10 |
| 4024 | SiVax Silicone Haematuria Catheter, 3-way, 24Fr. - Straight Round Tip | 10 |
| 4122 | SiVax Silicone Haematuria Catheter, 3-way, 22Fr. - Coude Round Tip | 10 |
| 4222 | SiVax Silicone Haematuria Catheter, 3-way, 22Fr. - Dufour Tip | 10 |
| 4322 | SiVax Silicone Haematuria Catheter, 3-way, 22Fr. - Couvelaire Tip | 10 |